
medisana®

## Infrared clinical thermometer NCT

### Unit and controls

- Backlight LCD display
- **START** button (starts reading and calls up stored data)
- **ON/OFF** button
- Battery compartment
- Lights
- Sensor
- Selector switch (measuring body / object temperature)

### Explanation of symbols

This instruction manual belongs to this device. It contains important information about starting up and operation. Read the instruction manual thoroughly. Non-observance of these instructions can result in serious injury or damage to the device.

**WARNING**  
These warning notes must be observed to prevent any injury to the user.

**CAUTION**  
These notes must be observed to prevent any damage to the device.

**NOTE**  
These notes give you useful additional information on the installation or operation.

**Classification:** Type BF applied part

**LOT** number

**Manufacturer**

**Date of manufacture**

## GB SAFETY INSTRUCTIONS

Read the instruction manual carefully before using this device, especially the safety instructions, and keep the instruction manual for future use. Should you give this device to another person, it is vital that you also pass on these instructions for use.

Only use the unit for the specific purpose described in the instruction manual. Any misuse will void the warranty. The unit is not designed for commercial use. Using this thermometer is no substitute for consulting a doctor. The thermometer and the patient must remain in a stable environment for at least 30 minutes before measuring the temperature. This device is not designed to be used by persons (including children) with limited physical, sensory or mental abilities, or by persons with insufficient experience and/or knowledge, unless under observation by a person responsible for their safety, or unless they have been instructed in the use of the device. Children must be supervised to ensure that they do not play with the device. The temperature in the forehead and temple area differs from the internal temperature, which is taken orally or rectally. Vasoconstriction, an effect which constricts the blood vessels and cools the skin, can occur during the early stages of a fever. In this case, the temperature measured by the NCT thermometer may be unusually low. If the measurement therefore does not match the patient's own perception or is unusually low, repeat the measurement every 15 minutes. As a reference, you can also measure the internal body temperature using a conventional oral or rectal thermometer. A control measurement using a conventional thermometer is recommended in the following cases:

- If the reading is surprisingly low
  - For new-born infants up to 100 days old
  - For children under three years of age who have a weakened immune system or who react unusually in the presence or absence of fever
- If you are using the thermometer for the first time or want to familiarise yourself with it
  - Babies' skin reacts very quickly to the ambient temperature. Therefore, do not take their temperature with the NCT thermometer during or after breastfeeding, because the skin temperature may then be lower than the internal body temperature.

If a fault occurs or the unit is damaged, do not attempt to repair the unit yourself, as this will void the warranty. Contact your dealer and only have repairs carried out by authorized service partners. Do not use the thermometer if it is damaged. Do not open up the device. It contains small parts which might be swallowed by children. Do not use the thermometer in damp environments. Do not immerse the thermometer in water or other liquids. It is not watertight. For cleaning, follow the instructions in **Cleaning and maintenance**. Do not shake or tap the thermometer. Do not drop it. Protect the device from direct sunlight, extremely high or low temperatures, dirt and dust. If the thermometer is not going to be used for a long period, remove the batteries.

**WARNING - The measurement results obtained by this thermometer do not represent a diagnosis and do not replace the consultation of a doctor, who should always be consulted in case of doubt.**

### Items supplied and packaging

Please check first of all that the device is complete and is not damaged in any way. In case of doubt, do not use it and contact your dealer or your service centre. The following parts are included:

- 1 **medisana** Infrared clinical thermometer NCT
- 2 Batteries (AAA type, LR03) 1.5V (already inserted!)
- 1 Storage box
- 1 Instruction manual

The packaging can be reused or recycled. Please dispose properly of any packaging material no longer required. If you notice any transport damage during unpacking, please contact your dealer without delay.

**WARNING - Please ensure that polythene packing is kept away from the reach of children! Risk of suffocation!**

### Taking a temperature reading - quickly and accurately!

The **medisana** NCT infrared clinical thermometer allows precise, contact-free measurement of body temperature (at a distance of approximately 5 cm). The temperature is taken above the eyebrow between the forehead and the temple. The value of the energy radiated there is measured and converted to a temperature value. The reading is displayed within 1 second. This enables you, for example, to take babies' temperature without disturbing their sleep.

### Additional function

Apart from measuring body temperature, the NCT infrared thermometer from **medisana** can also be used to measure the temperature of objects (e.g. fluids, surfaces, ambient temperature). The temperature range for this function is between 0° and 100° Celsius (32°F-212°F).

## BATTERY SAFETY INFORMATION

Do not disassemble batteries!  
Clean the battery and device contacts if necessary before putting in the batteries!  
Remove discharged batteries from the device immediately!  
Increased risk of leakage, avoid contact with skin, eyes and mucous membranes! If battery acid comes in contact with any of these parts, rinse the affected area with copious amounts of fresh water and seek medical attention immediately!  
If a battery has been swallowed seek medical attention immediately!  
Replace all of the batteries simultaneously!  
Only replace with batteries of the same type, never use different types of batteries together or used batteries with new ones!  
Insert the batteries correctly, observing the polarity!  
Keep the battery compartment well sealed!  
Remove the batteries from the device if it is not going to be used for an extended period!  
Keep batteries out of children's reach!  
Do not attempt to recharge these batteries! There is a danger of explosion!  
Do not short circuit! There is a danger of explosion!  
Do not throw into a fire! There is a danger of explosion!  
Keep unused batteries in their packaging away from metal objects in order to prevent short circuiting!  
Do not throw used batteries into the household refuse; put them in a hazardous waste container or take them to a battery collection point, at the shop where they were purchased!

### Removing / inserting batteries

The **medisana** NCT infrared thermometer comes with the batteries installed. However, if the battery symbol appears in the LC display after a certain period of use, you must replace them with two new 1.5 V batteries (LR03, AAA).

#### Removal:

Slide the lid of the battery compartment cover downwards and take the cover off. Take out the batteries.

#### Insertion:

Put in two new 1.5 V batteries (LR03, AAA). Make sure the batteries are the right way around (see the diagram in the battery compartment). Put the battery compartment cover back on and move it upwards until it locks into position.

### Before use

Acclimatise the thermometer to the temperature of the room: For an accurate reading of a patient's temperature, their body must first acclimatise to the temperature of the room where the measurement will be made. The same applies to the NCT infrared thermometer. When taken somewhere else, it automatically adjusts itself to the new environment. The thermometer takes a certain amount of time to equalise its own temperature to that of the room. This depends on the difference between the two temperatures and can take more than 30 minutes. Precise readings can only be taken once this process is completed.

### Change the unit of measurement

You may measure the body temperature in °C (Celsius) or °F (Fahrenheit). This is how to do it:

The unit must be switched off.

Press and hold the **START**-button for at least 5 seconds, until the current selected unit of measurement is displayed. Another short press on the **START**-button changes the unit of measurement. When the desired unit is displayed, wait until all symbols are displayed. Now the thermometer is ready for the measurement in the new selected unit of measurement.

Additionally all saved measurements are adjusted automatically according to the newly selected unit of measurement.

### Measuring body temperature

**NOTES, for the most accurate readings:**

The thermometer and the patient must remain in a stable environment for at least 30 minutes before measuring the temperature.

Do not take a baby's temperature during or directly after breastfeeding.

Do not use the thermometer in damp environments.

Patients should not drink or eat directly before or during a measurement, and they should avoid any physical exertion.

Clean the area where you are taking the temperature and remove any dirt, hairs or sweat before moving the sensor to the measuring position.

Do not remove the thermometer from the measuring position until it has beeped to signal the end of the reading.

If the sensor becomes dirty during or after a reading, clean it with a cotton bud dipped in cleaning alcohol before putting the thermometer away.

Always take the temperature in the same place, because the readings will vary in different environments.

A control measurement using a conventional thermometer is recommended in the following cases:

- If the reading is surprisingly low,
- For new-born infants up to 100 days old,
- For children under three years of age who have a weakened immune system or who react unusually in the presence or absence of fever,
- If you are using the thermometer for the first time or want to familiarise yourself with it.

The thermometer has been clinically tested, and its safety and accuracy inspected and confirmed, provided it is used in accordance with these instructions.

Move the selector switch to the top position to select body temperature measurement mode .

Press the **ON/OFF** button . The LC display is activated. The thermometer performs a self-test, where all the display symbols are shown in succession.

Finally, the last registered reading is displayed along with the M symbol (M = memory) for 2 seconds.

When you hear a short beep and the symbol for the measuring unit °C is shown, the thermometer is ready for use.

To begin the measurement, point the sensor at the middle of the patient's forehead at a distance of no more than 5 cm. If the patient's forehead has hair, sweat or dirt on it, remove this first so that the reading is as accurate as possible.

Press and hold the **START** button . The blue light switches on to illuminate the area where you are taking the measurement.

Now move the thermometer smoothly over the patient's eye-brow from the forehead to the temple to measure the maximum temperature. While the sensor scans this area, you will hear a long beep to indicate that the maximum temperature has been recorded.

As soon as this happens, release the **START** button .

- If the reading is below the fever level, you will hear a long beep and the backlight of the display will be green.

- If the reading is at or above the fever level (37.8°C resp. 100°F), you will hear 10 short beeps, and the backlight of the display will be red to indicate that the patient may have a fever.

When the °C symbol flashes again after 2 seconds, the thermometer is ready for another reading.

### Measuring the temperature of an object

Move the selector switch to the bottom position to select object temperature measurement mode .

Press the **ON/OFF** button . The LC display is activated. The thermometer performs a self-test, where all the display symbols are shown in succession.

Finally, the last registered reading is displayed along with the M symbol (M = memory) for 2 seconds.

When you hear a short beep and the symbol for the measuring unit °C flashes, the thermometer is ready to take the measurement.

To begin the measurement, point the sensor at the middle of the object at a distance of no more than 5 cm. If it has water, dust or dirt on it, remove this first so that the reading is as accurate as possible.

Press the **START** button . The blue light switches on and the temperature is measured. The measurement is accompanied by a series of beeps.

Release the **START** button after a few seconds. A long beep indicates the end of the measurement. Remove the sensor from the measuring position. The display shows the measured temperature.

### Recalling stored data

The **medisana** NCT infrared thermometer can store up to 30 results. You can call up your stored data with a short press of the **START** button , without having previously switched on the thermometer with the **ON/OFF** button . The symbol M will flash on the display next to "----". After a further short press of the **START** button , memory slot number 1 will appear in the unit display, followed directly by the last stored temperature reading. The symbol M will flash below the unit of measurement. With a further press of the **START** button , you can call up the result stored in memory slot "2". A total of 30 memory slots are available, and these will be displayed in order if you continue to press the **START** button . If a memory slot is empty, the symbol "----" will appear next to it in the display.

### Error messages

Display	Description
	The measured temperature is too high. H appears in the display if the reading is above 40°C / 108°F (body mode) or 100°C / 212°F (object mode).
	The ambient temperature is too low. L appears in the display if the reading is below 34°C / 93.2°F (body mode) or 0°C / 32°F (object mode).
	The ambient temperature is too high. The ambient temperature is above 40°C / 104°F (body and object mode).
	The ambient temperature is too low. The ambient temperature is below: 15°C / 59°F (body mode) or 5°C / 41°F (object mode).
	Error indication in the display. A malfunction has occurred. You hear 3 short beeps.
	Blank display. Check the position of the battery.
	Flat battery indicated. Insert a new battery.

### Cleaning and maintenance

- The sensor is the most delicate part of the **medisana** NCT infrared thermometer. Protect it from dirt and damage. Clean the sensor with a cotton bud dipped in warm water or 70% isopropyl alcohol.
- Do not immerse the thermometer in water or any other liquid. It is not watertight.
- To clean the thermometer, use a soft cloth slightly dampened with mild soapy water or 70% isopropyl alcohol. Do not let water get into the unit. Only use the unit again when completely dry.
- Never use aggressive cleaning agents, thinners, benzine or tough brushes.
- Take care not to scratch the lens of the sensor, the display or the housing.
- Do not use or store the unit in an environment where the ambient temperature or humidity is too high or too low (see operating conditions in the "Technical specifications" section), in direct sunshine, in the vicinity of static electricity or in a dusty environment.
- Remove the batteries if you do not wish to use the unit for an extended period of time. Otherwise there is a danger that the batteries may leak.
- If the thermometer is used as described in these instructions, no regular maintenance or calibration is required.

### Disposal

This product must not be disposed of together with domestic waste. All users are obliged to hand in all electrical or electronic devices, regardless of whether or not they contain toxic substances, at a Z municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner. Please remove the batteries before disposing of the device/unit. Do not dispose of old batteries with your household waste, but at a battery collection station at a recycling site or in a shop. Consult your municipal authority or your dealer for information about disposal.

### Technical specifications

Name and model:	<b>medisana</b> Infrared clinical thermometer NCT
Display system:	digital display (LCD)
Power supply:	3 V <sub>z</sub> , 2 batteries 1.5 V (type LR03, AAA)
Measuring range:	Body mode: 34°C - 42.2°C (93.2°F - 108°F) Object mode: 0°C - 100°C (32°F - 212°F) in the region of 34°C - 35.9°C: ±0.3°C (93.2°F - 96.6°F): ±0.54°F 36°C - 39°C: ±0.2°C (96.8°F - 102.2°F): ±0.36°F 39.1°C - 42.2°C: ±0.3°C (102.4°F - 108°F): ±0.54°F (under operating conditions) 0.18°C (± 0.324°F)
Accuracy of measurement: (laboratory readings)	approx. 1 second up to 30 readings for temperatures > 37.8°C (100°F) approx. 3 minutes after the last measurement Body mode 15°C - 40°C (59°F - 104°F) Object mode 5°C - 40°C (41°F - 104°F)
Clinical repeatability:	approx. 147 x 38 x 21 mm
Duration of measurement:	approx. 48 g including batteries
Memory capacity:	48614
Fever alarm:	40 15588 48614 2
Autom. switch-off:	
Operating temperature:	
Storage/transport temperature:	
Dimensions:	
Weight:	
Article No.:	
EAN Code:	

The unit is designed for home use only. For commercial use, a metrological inspection must be carried out once every year. This calibration will be charged for and can be carried out by an appropriate authority or an authorised service centre - in accordance with the regulations for users of medical equipment.

### Electromagnetic compatibility:

The device complies with the EN 60601-1-2 standard for electro-magnetic compatibility. You can find details on this measurement data on the **medisana** homepage [www.medisana.com](http://www.medisana.com).

In accordance with our policy of continual product improvement, we reserve the right to make technical and optical changes without notice.

### Warranty and repair terms

Please contact your dealer or the service centre in case of a claim under the warranty. If you have to return the unit, please enclose a copy of your receipt and state what the defect is.


The following warranty terms apply:

1. The warranty period for **medisana** products is three years from date of purchase. In case of a warranty claim, the date of purchase has to be proven by means of the sales receipt or invoice.
2. Defects in material or workmanship will be removed free of charge within the warranty period.
3. Repairs under warranty do not extend the warranty period either for the unit or for the replacement parts.
4. The following is excluded under the warranty:
  - a. All damage which has arisen due to improper treatment, e.g. nonobservance of the user instructions.
  - b. All damage which is due to repairs or tampering by the customer or unauthorised third parties.
  - c. Damage which has arisen during transport from the manufacturer to the consumer or during transport to the service centre.
  - d. Accessories which are subject to normal wear and tear.
5. Liability for direct or indirect consequential losses caused by the unit are excluded even if the damage to the unit is accepted as a warranty claim.

**Microlife Corporation**  
9F, 431, RuiGuang Road, NeiHu  
Taipei 11492, Taiwan, China

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**Items supplied and packaging**  
Please check first of all that the device is complete and is not damaged in any way. In case of doubt, do not use it and contact your dealer or your service centre. The following parts are included:

- 1 **medisana** Infrared clinical thermometer **NCT**
- 2 Batteries (AAA type, LR03) 1.5V (already inserted!)
- 1 Storage box
- 1 Instruction manual

The packaging can be reused or recycled. Please dispose properly of any packaging material no longer required. If you notice any transport damage during unpacking, please contact your dealer without delay.

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