
medisana®

Infrared Ear & Forehead Thermometer

TM 300

GB Unit and controls

- 1 Sensor
- 2 SCAN-button
- 3 Display
- 4 Battery compartment (at the back side)
- 5 MODE-button
- 6 Sensor cover
- 7 Ambient measurement
- 8 Forehead measurement
- 9 In-ear measurement
- 10 Time
- 11 Fever
- 12 Normal temperature
- 13 Measurement unit (°C or °F)
- 14 Memory display
- 15 Battery symbol

Explanation of symbols

This instruction manual belongs to this device. It contains important information about starting up and operation. Read the instruction manual thoroughly. Non-observance of these instructions can result in serious injury or damage to the device.

WARNING
These warning notes must be observed to prevent any injury to the user.

CAUTION
These notes must be observed to prevent any damage to the device.

NOTE
These notes give you useful additional information on the installation or operation.

Classification: Type BF applied part

Temperature range

Keep dry

IP22 Protection classification against solid objects and water

LOT LOT number

Manufacturer

Date of manufacture

GB Infrared Ear & Forehead Thermometer

Manual - Please read carefully!

Safety Information

IMPORTANT INFORMATION! RETAIN FOR FUTURE USE!

Read the instruction manual carefully before using this device, especially the safety instructions, and keep the instruction manual for future use. Should you give this device to another person, it is vital that you also pass on these instructions for use.

INTENDED USE:

Thank you for purchasing our Infrared Ear & Forehead Thermometer TM 300. This Infrared Thermometer is a non-sterile, reusable, handheld device. It can be used by consumers in home care environment and doctors in clinic as reference. It is intended for measuring human body temperature of three months and above of people by detecting infrared heat from the center of the forehead and in the ear.

- Only use the unit for the specific purpose described in the instruction manual. Any misuse will void the warranty.
- The unit is not designed for commercial use.
- Using this thermometer is no substitute for consulting a doctor.
- The thermometer and the patient must remain in a stable environment for at least 30 minutes before measuring the temperature.
- This device is not designed to be used by persons (including children) with limited physical, sensory or mental abilities, or by persons with insufficient experience and/or knowledge, unless under observation by a person responsible for their safety, or unless they have been instructed in the use of the device.
- Children must be supervised to ensure that they do not play with the device.
- The temperature in the forehead and temple area differs from the internal temperature, which is taken orally or rectally. If the measurement therefore does not match the patient's own perception or is unusually low, repeat the measurement every 15 minutes. As a reference, you can also measure the internal body temperature using a conventional oral or rectal thermometer.
- A control measurement using a conventional thermometer is recommended in the following cases:
 1. If the reading is surprisingly low
 2. For new-born infants up to 100 days old
 3. For children under three years of age who have a weakened immune system or who react unusually in the presence or absence of fever
 4. If you are using the thermometer for the first time or want to familiarise yourself with it
- If a fault occurs or the unit is damaged, do not attempt to repair the unit yourself, as this will void the warranty. Contact your dealer and only have repairs carried out by authorized service partners.
- Do not modify and / or open the device. It contains small parts which might be swallowed by children.
- Do not use the thermometer if it is damaged.
- Do not use the thermometer in damp environments. Do not immerse the thermometer in water or other liquids. **It is not watertight.**
- This device comprises sensitive components and must be treated with caution. Observe the storage and operating conditions described in the «Technical Specifications» section.
- Do not use this device close to strong electromagnetic fields such as mobile telephones or radio installations.
- Keep a minimum distance of 3.3 m from such devices when using this device.
- Do not shake or tap the thermometer. Do not drop it.
- Protect the device from direct sunlight, extremely high or low temperatures, dirt and dust.
- If the thermometer is not going to be used for a long period, remove the battery.

WARNING - The measurement results obtained by this thermometer do not represent a diagnosis and do not replace the consultation of a doctor, who should always be consulted in case of doubt.

Items supplied and packaging

Please check first of all that the device is complete and is not damaged in any way. In case of doubt, do not use it and contact your dealer or your service centre. The following parts are included:

- 1 **medisana** Infrared Ear & Forehead Thermometer **TM 300**
- 1 Storage bag
- 2 Batteries (AAA type, LR03) 1.5V
- 1 Instruction manual

The packaging can be reused or recycled. Please dispose properly of any packaging material no longer required. If you notice any transport damage during unpacking, please contact your dealer without delay.

WARNING - Please ensure that polythene packing is kept away from the reach of children! Risk of suffocation!

Taking a temperature reading - quickly and accurately!

The **medisana** Infrared Ear & Forehead Thermometer allows precise measurement of body temperature. The temperature is taken above the eyebrow between the forehead and the temple.

The value of the energy radiated there is measured and converted to a temperature value.

Additional function: Apart from measuring body temperature, the Infrared multifunctional thermometer can also be used to measure the temperature of objects (e.g. fluids, surfaces, ambient temperature). The temperature range for this function is between 0-100°C (32-212°F).

Operation - Getting started

When the battery symbol appears in the display after a longer time of usage, the batteries should be replaced soonest. When the battery symbol appears or nothing at all appears in the display , the batteries need to be replaced immediately.

Removal

Remove the lid of the battery compartment cover and take out the batteries.

Insertion

Put in two new 1.5V batteries of type LR03/AAA. Make sure the batteries are the right way around (see the diagram in the battery compartment). Put the battery compartment lid back.

BATTERY SAFETY INFORMATION

- Do not disassemble batteries!
- Clean the battery and device contacts if necessary before putting in the batteries!
- Remove discharged batteries from the device immediately!
- Increased risk of leakage, avoid contact with skin, eyes and mucous membranes! If battery acid comes in contact with any of these parts, rinse the affected area with copious amounts of fresh water and seek medical attention immediately!
- If a battery has been swallowed seek medical attention immediately!
- Insert the battery correctly, observing the polarity!
- Keep the battery compartment well sealed!
- Remove the battery from the device if it is not going to be used for an extended period!
- Keep batteries out of childrens reach!
- Do not attempt to recharge batteries! There is a danger of explosion!
- Do not short circuit! There is a danger of explosion!
- Do not throw into a fire! There is a danger of explosion!
- Keep unused batteries in their packaging away from metal objects in order to prevent short circuiting!
- Do not throw used batteries into the household refuse; put them in a hazardous waste container or take them to a battery collection point, at the shop where they were purchased!

Settings

After inserting a battery, all elements of the display are shown for approx. 1 second. Afterwards the device switches to the setting mode (the display shows the flashing hour figure "12"). Use the SCAN-button to adjust the hour and press the MODE-button to confirm. Now adjust minute in the same way and press the MODE-button to confirm the setting. Afterwards, the device is ready for measurement.

NOTES

- The device will automatically switch to time display if there is no input for approx. 25 seconds.

Before use

Acclimatise the thermometer to the temperature of the room: For an accurate reading of a patient's temperature, their body must first acclimatise to the temperature of the room where the measurement will be made. The same applies to the Infrared Ear & Forehead Thermometer TM 300. When taken somewhere else, it automatically adjusts itself to the new environment. The thermometer takes a certain amount of time to equalise its own temperature to that of the room. This depends on the difference between the two temperatures and can take more than 30 minutes. Precise readings can only be taken once this process is completed.

Measuring forehead temperature

NOTES

- For the most accurate readings:
 - Do not take a baby's temperature during or directly after breastfeeding.
 - Patients should not drink or eat directly before or during a measurement, and they should avoid any physical exertion.
 - Clean the area where you are taking the temperature and remove any dirt, hairs or sweat before moving the sensor to the measuring position.
 - Do not remove the thermometer from the measuring position until it has beeped to signal the end of the reading.
 - If the sensor becomes dirty during or after a reading, clean it with a cotton bud dipped in cleaning alcohol before putting the thermometer away.
 - Always take the temperature in the same place, because the readings will vary in different environments.

The thermometer's safety and accuracy has been inspected and confirmed, provided it is used in accordance with these instructions.

1. Attach the sensor cover to the sensor before measurement. Then press the SCAN-button to activate the thermometer and the last measurement is displayed. After approx. 3 seconds, the measurement unit flashes and the device is ready for measurement - this will be indicated by two short beeps.
2. Press the MODE-button and switch to Forehead measurement .
3. To begin the measurement, put the sensor at the middle of the patient's forehead. If the patient's forehead has hair, sweat or dirt on it, remove this first so that the reading is as accurate as possible.
4. Press the SCAN-button and release. After approx. 1 second, a short beep can be heard to indicate that the measurement is finished and the measured temperature will be shown.
5. If the result is 37.5-37.9°C, display will show a smiling face with green color backlight, which means normal temperature. If the result is $\geq 37.5^\circ\text{C}$, a smiling face and orange color backlight will be shown, that means fever is detected. Backlight will turn to red color with sad face if the result is $\geq 38^\circ\text{C}$.

If temperature $\leq 32^\circ\text{C}$ (89.6°F), «Lo» will display. If temperature $\geq 42.9^\circ\text{C}$ (109.2), «Hi» will display.

After approx. 5 seconds, the measurement unit flashes and the device is ready for another measurement. The device will automatically switch off if there is no input for approx. 25 seconds and will automatically switched to time display.

Measuring in-ear temperature

1. To measure the temperature in your ear, first remove the sensor cover . Then press the SCAN-button to activate the thermometer and the last measurement is displayed. After approx. 3 seconds, the measurement unit flashes and the device is ready for measurement - this will be indicated by two short beeps.
2. Stretch your ear canal by pulling your ear backwards and upwards and carefully insert the sensor into the ear. For children under the age of 1 year, the sensor should be inserted into the ear canal, while the child is lying flat (e.g. on a bed) with the head slightly turned, so that the ear points upward.
3. Press the SCAN-button and release. After approx. 1 second a long beep can be heard to indicate that the measurement is finished and the measured temperature will be shown.
4. If the result is $\leq 37.4^\circ\text{C}$, display will show a smiling face with green color backlight, which means normal temperature. If the result is 37.5-37.9°C, a smiling face and orange color backlight will be shown, that means fever is detected. Backlight will turn to red color with sad face if the result is $\geq 38^\circ\text{C}$.

If temperature $\leq 32^\circ\text{C}$ (89.6°F), «Lo» will display. If temperature $\geq 42.9^\circ\text{C}$ (109.2), «Hi» will display.

After approx. 5 seconds, the measurement unit flashes and the device is ready for another measurement. The device will automatically switch off if there is no input for approx. 25 seconds and will automatically switched to time display.

Measuring object temperature

1. Attach the sensor cover to the sensor before measurement. Then press the SCAN-button to activate the thermometer and the last measurement is displayed. After approx. 3 seconds, the measurement unit flashes and the device is ready for measurement - this will be indicated by two short beeps.
2. Press the MODE-button for twice and switch to Object measurement (no symbol will be indicated).
3. To begin the measurement, put the sensor at the middle of the object. Press the SCAN-button and release. After approx. 1 second, a short beep can be heard to indicate that the measurement is finished and the measured temperature will be shown.

After approx. 5 seconds, the measurement unit flashes and the device is ready for another measurement. The device will automatically switch off if there is no input for approx. 25 seconds and will automatically switched to time display.

Change the measurement unit

In clock display, long-press the SCAN-button until you see the measurement unit (°C or °F) selected. Then press SCAN-button to start measurement.

Recalling stored data

The thermometer can store up to 20 results. You can call up your stored data with press and release the SCAN-button . And then, long press the MODE-button until «00» shows in display. If measurement values are stored in memory, first the memory slot number «00» appears with the stored temperature. By further pressing the MODE-button , you can call up more stored values (memory slot «02» etc.). Totally there are 20 memory slots available, which may be recalled by pressing the MODE-button repeatedly.

Error messages

Screen display	Reasons	Solutions
HH	When the environment temperature is higher than normal operation measurement range, the LCD screen will show "HH" indication.	When you operate this product, the environment temperature must be 15-40°C. Please place this product and keep steady in current environment for 30 minutes before next measurement.
LL	When the environment temperature is lower than normal operation measurement range, the LCD screen will show "LL" indication.	
Hi	When the target object temperature is higher than the limit of the device.	
Lo	When the target object temperature is higher than the limit of the device.	

Troubleshooting

Phenomenons	Reasons	Solutions
When the power is turned on, the screen can not be displayed.	The battery is exhausted.	Replace the new battery.
	The battery polarity is wrong.	The battery polarity is the same as the battery case.
The measurement temperature is low.	The measurement position is not correct.	Measure the temperature correctly according to the instructions.
	There are dirt stuff in the sensor or ear canal.	Please clear the dirt before the measurement.
Big temperature fluctuations with continuous measurement.	The measurement interval is too small.	The interval for each measurement should be above 10 seconds.

Cleaning and maintenance

- The sensor is the most delicate part of the device. Protect it from dirt and damage. Clean the sensor with a cotton bud dipped in 70% alcohol.
- Do not immerse the thermometer in water or any other liquid. **It is not watertight.**
- To clean the thermometer, use a soft cloth slightly dampened with 70% alcohol. Do not let water get into the unit. Only use the unit again when completely dry.
- Never use aggressive cleaning agents, thinners, benzine or tough brushes.
- Take care not to scratch the lens of the sensor, the display or the housing.
- Do not use or store the unit in an environment where the ambient temperature or humidity is too high or too low (see operating conditions in the "Technical specifications" section), in direct sunshine, in the vicinity of static electricity or in a dusty environment.
- Remove the batteries if you do not wish to use the unit for an extended period of time. Otherwise there is a danger that the batteries may leak.
- If the thermometer is used as described in these instructions, no regular maintenance or calibration is required.

Disposal

This product must not be disposed of together with domestic waste. All users are obliged to hand in all electrical or electronic devices, regardless of whether or not they contain toxic substances, at a municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner. Please remove the battery before disposing of the device/unit. Do not dispose of old batteries with your household waste, but at a battery collection station at a recycling site or in a shop. Consult your municipal authority or your supplier for information about disposal.

Technical specifications

Name and model:	medisana Infrared Ear & Forehead Thermometer TM 300
Display system:	digital display (LCD)
Power supply:	2 x batteries (AAA type, LR03) 1.5V
Expected service life:	5 years or 12,000 measurements
Measuring range:	In-ear and forehead mode 32-42.9°C (89.6-109.2°F); Object and ambient mode 0-100°C (32-212°F)
Accuracy of measurement (laboratory readings):	In-ear and forehead mode 32-42.9°C (89.6-109.2°F); Object and ambient mode 0-100°C (32-212°F) In-ear and forehead mode: $\pm 0.2^\circ\text{C}$, 35-42°C ($\pm 0.4^\circ\text{F}$, 95-107.6°F); 0.1°C (0.1°F)
Resolution:	up to 20 readings
Memory capacity:	5-40°C (41-104°F); rel. humidity 15-95% RH
Operating temperature:	-25°C to +55°C (-13°F to +131°F), rel. humidity 15-95%
Storage/transport temperature:	approx. 145 x 37 x 35 mm
Dimensions:	approx. 53 g
Weight:	48604
Article No.:	48604
EAN Code:	40155 88 48604 3

The CE-mark of this device refers to the EC-directive 93/42/EEC MDD. The unit is designed for home use only. For commercial use, a metrological inspection must be carried out once every year. This calibration will be charged for and can be carried out by an appropriate authority or an authorised service centre - in accordance with the regulations for users of medical equipment.

Electromagnetic compatibility: The device complies with the EN 60601-1-2 standard for electromagnetic compatibility.

In accordance with our policy of continual product improvement, we reserve the right to make technical and optical changes without notice.

Warranty and repair terms

Please contact your dealer or the service centre in case of a claim under the warranty. If you have to return the unit, please enclose a copy of your receipt and state what the defect is. The following warranty terms apply.

1. The warranty period for **medisana** products is three years from date of purchase. In case of a warranty claim, the date of purchase has to be proven by means of the sales receipt or invoice.
2. Defects in material or workmanship will be removed free of charge within the warranty period.
3. Repairs under warranty do not extend the warranty period either for the unit or for the replacement parts.
4. The following is excluded under the warranty:
 - a. All damage which has arisen due to improper treatment, e.g. nonobservance of the user instructions.
 - b. All damage which is due to repairs or tampering by the customer or unauthorised third parties.
 - c. Damage which has arisen during transport from the manufacturer to the consumer or during transport to the service centre.
 - d. Accessories which are subject to normal wear and tear.
5. Liability for direct or indirect consequential losses caused by the unit are excluded even if the damage to the unit is accepted as a warranty claim.

Kangfu Medical Equipment Factory
No. 380 Ningkang East Road, Yueqing City, P.R.C.

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EC REP Luxus Lebenswelt GmbH