
medisana®

Non-contact Infrared Thermometer FTE Instruction manual

Explanation of symbols

IMPORTANT
Non-observance of these instructions can result in serious injury or damage to the device.

WARNING
These warning notes must be observed to prevent any injury to the user.

CAUTION
These notes must be observed to prevent any damage to the device.

NOTE
These notes give you useful additional information on the installation and operation.

Device classification: type BF

Lot number

Manufacturer

Date of manufacture

48610_FTE_manual_GB_202105_L5

CE 0297

GB Safety Instructions

Important Information! Retain for future use! Read the instruction manual carefully before using this device, especially the safety instructions, and keep the instruction manual for future use. Should you give this device to another person, it is vital that you also pass on these instructions for use.

Intended use:

The **medisana** infrared clinical thermometer **FTE** is suitable to measure the body temperature of humans. The measurement is carried out contact-free on the forehead. Additionally the thermometer is suitable to measure the temperature of objects in the range of 0°C to 100°C.

- Only use the unit for the specific purpose described in the instruction manual. Any misuse will void the warranty.
- The unit is not designed for commercial use.
- Using this thermometer is no substitute for consulting a doctor. Consult your doctor in any case, before you start a medicamentous treatment.
- The thermometer and the patient must remain in a stable environment for at least 30 minutes before measuring the temperature. Do not take a measurement within 30 minutes after physical exercise or after having a bath.
- This device is not designed to be used by persons (including children) with limited physical, sensory or mental abilities, or by persons with insufficient experience and/or knowledge, unless under observation by a person responsible for their safety, or unless they have been instructed in the use of the device.
- Children must be supervised to ensure that they do not play with the device.
- The temperature in the forehead and temple area differs from the internal temperature, which is taken orally or rectally. Vasoconstriction, an effect which constricts the blood vessels and cools the skin, can occur during the early stages of a fever. In this case, the temperature measured by the **infrared clinical thermometer FTE** may be unusually low. If the measurement therefore does not match the patient's own perception or is unusually low, repeat the measurement every 15 minutes. As a reference, you can also measure the internal body temperature using a conventional oral or rectal thermometer.
- It is recommended to perform at least 3 measurements with at least 1 minute interval. If the readings are different, accept the highest value as reference value.
- A control measurement using a conventional thermometer is recommended in the following cases:
 - If the reading is surprisingly low
 - For newborn infants up to 100 days old
 - For children under three years of age who have a weakened immune system or who react unusually in the presence or absence of fever
 - If you are using the thermometer for the first time or want to familiarise yourself with it

- Babies' skin reacts very quickly to the ambient temperature. Therefore, do not take their temperature with the **infrared clinical thermometer FTE** during or after breastfeeding, because the skin temperature may then be lower than the internal body temperature.
- There are no absolute body temperature standards. Keep record of your personal readings to allow your doctor to judge if a medical treatment is necessary or not.
- Do not change the construction of the device and do not attempt to repair the unit yourself if a fault occurs or the unit is damaged, as this will void the warranty. Contact your dealer and only have repairs carried out by authorized service partners.
- Do not use the thermometer if it is damaged.
- Do not shake or tap the thermometer. Do not drop it.
- Do not use the thermometer in damp environments.
- Do not immerse the thermometer in water or other liquids. **It is not watertight.**
- The sensor should be cleaned after each use, see the instructions in **Cleaning and maintenance**.
- Avoid direct finger contact with the sensor.
- Do not shake or tap the thermometer. Do not drop it.
- Protect the device from direct sunlight, extremely high or low temperatures, dirt and dust.
- If the thermometer is not going to be used for a long period, remove the batteries.

Unit and controls

- Illuminated LCD display
- Signal LED (fever alarm)
- MODE-button
- MEM-button
- SCAN-button
- Sensor
- Battery compartment cover
- Symbol for running measurement
- Night mode

- Measurement of body temperature (Body mode)
- Measurement of object temperature (Object mode)
- Temperature display
- °C Celsius
- °F Fahrenheit
- Battery exchange symbol
- Fever alarm symbol ($\geq 38.0^{\circ}\text{C}/100.4^{\circ}\text{F}$)
- Symbol for normal temperature ($< 38.0^{\circ}\text{C}/100.4^{\circ}\text{F}$)
- Memory display

Items supplied and packaging

Please check first of all that the device is complete and is not damaged in any way. If in doubt, do not use it. Send it to a service point.

The following parts are included:

- 1 **medisana** Infrared clinical thermometer **FTE**
- 2 Batteries (AAA type, LR03) 1.5V
- 1 Instruction manual

The packaging can be reused or recycled. Please dispose properly of any packaging material no longer required. If you notice any transport damage during unpacking, please contact your dealer without delay.

WARNING
Please ensure that polythene packing is kept away from the reach of children! Risk of suffocation!

Taking a temperature reading - quickly and accurately!

The **medisana** Infrared clinical thermometer allows precise, contact-free (approx. 0.5 to 3 cm distance) measurement of body temperature. The temperature is taken above the eyebrow between the forehead and the temple. The value of the energy radiated there is measured and converted to a temperature value. The reading is displayed within a few seconds. This enables you, for example, to take babies' temperature without disturbing their sleep.

Additional function

Apart from measuring body temperature, the **medisana** IR-clinical thermometer can also be used to measure the temperature of objects (e.g. fluids, surfaces, ambient temperature). The temperature range for this function is between 0°C and 100°C (32°F - 212°F).

Getting started: Removing / Inserting batteries


For the use resp. if the battery exchange symbol appears in the display after a certain period of use, you must insert new batteries. If the battery exchange symbol appears together with "Lo" in the display, the batteries need to be replaced before taking another measurement. Use two new 1.5 V batteries (LR03, AAA).

Removal: Take off the battery compartment cover

by sliding it downwards as shown. Take out the batteries.

Insertion: Put in two new 1.5 V batteries (LR03, AAA).

Make sure the batteries are the right way around (see the diagram in the battery compartment). Put the battery compartment cover back on.

WARNING BATTERY SAFETY INFORMATION

- Do not disassemble batteries!
- Clean the battery and device contacts if necessary before putting in the batteries!
- Remove discharged batteries from the device immediately!
- Increased risk of leakage, avoid contact with skin, eyes and mucous membranes! If battery acid comes in contact with any of these parts, rinse the affected area with copious amounts of fresh water and seek medical attention immediately!
- If a battery has been swallowed seek medical attention immediately!
- Replace all of the batteries simultaneously!
- Only replace with batteries of the same type, never use different types of batteries together or used batteries with new ones!
- Insert the batteries correctly, observing the polarity!
- Keep the battery compartment well sealed!
- Remove the batteries from the device if it is not going to be used for an extended period!
- Keep batteries out of children's reach!
- Do not attempt to recharge these batteries! There is a danger of explosion!
- Do not short circuit! There is a danger of explosion!
- Do not throw into a fire! There is a danger of explosion!
- Keep unused batteries in their packaging away from metal objects in order to prevent short circuiting!
- Do not throw used batteries into the household refuse; put them in a hazardous waste container or take them to a battery collection point, at the shop where they were purchased!

Before use

NOTE, for the most accurate readings:

- The thermometer and the patient must remain in a stable environment for at least 30 minutes before measuring the temperature.
- Do not take a baby's temperature during or directly after breastfeeding.
- Do not use the thermometer in damp environments.
- Patients should not drink or eat directly before or during a measurement, and they should avoid any physical exertion.
- Clean the area where you are taking the temperature and remove any dirt, hairs or sweat before moving the sensor to the measuring position.
- Do not remove the thermometer from the measuring position until it has beeped to signal the end of the reading.
- If the sensor becomes dirty during or after a reading, clean it with a cotton bud dipped in cleaning alcohol before putting the thermometer away.
- Always take the temperature in the same place, because the readings will vary in different environments.
- The thermometer has been clinically tested, and its safety and accuracy inspected and confirmed, provided it is used in accordance with these instructions.

Change the unit of measure


You can measure the temperature in ° Celsius or in ° Fahrenheit. Please proceed as follows to select the unit of measure:

- The device has to be switched off.
- Press the **SCAN**-button to switch on the device. Shortly all display elements are shown (initializing the display).
- Press the **MODE**-button and the **MEM**-button simultaneously and hold for approx. 3 seconds. This will change the unit of measure to either °C or °F. When the selected unit appears in the display, wait until the device shortly beeps to confirm the selection. After that the unit is ready for the next measurement in the selected unit of measure.

Select the measurement mode

You can select one out of 4 measurement modes. Please proceed as follows to select the measurement mode:

- When the unit switched on, press the **MODE**-button to select one of the below four measurement modes:

- Press the **MODE**-button until the desired measurement mode is displayed.

NOTE

The night mode deactivates the key tones to avoid unnecessary noise.

Measuring the body temperature

- Select a measurement mode (please refer to above chapter): Body mode or Body mode with night mode.
- To start the measurement, point the sensor in the middle of the patient's forehead at a distance between 0.5 and max. 3 cm. If the patient's forehead has hair, sweat or dirt on it, remove this first so that the reading is as accurate as possible.
- Press and hold the **SCAN**-button and move the thermometer smoothly over the patient's eyebrow from the forehead to the temple to measure the maximum temperature.
- Release the **SCAN**-button. You will hear a short beep.
- The measurement result will be shown on the display. After approx. 30 seconds without any further use the thermometer will switch off automatically to save battery power.

NOTE

Is the measured temperature below 38.0°C (100.4°F), a smiley face and a beep will be emitted. Is the measured temperature equal to or higher than 38.0°C (100.4°F), a frowning face and three beeps will be emitted and additionally the signal LED will light up in red.

Measuring the temperature of an object

- Select a measurement mode (please refer to above chapter): Object mode or object mode with night mode.
- To start the measurement, point the sensor on the object with a distance of 0.5 to max. 3 cm.
- Press and hold the **SCAN**-button for a few seconds.
- Release the **SCAN**-button. You will hear a short beep.
- The measurement result will be shown on the display. After approx. 30 seconds without any further use the thermometer will switch off automatically to save battery power.

Recalling stored data

The **medisana** Infrared clinical thermometer can store up to 30 results. You can call up your stored data with a short press on the **MEM**-button, while the appliance is switched on. The last stored temperature reading will be shown in the display, together with the memory slot number (slot number 1) and the symbol for the measurement mode. With each further press on the **MEM**-button you can call up the previous results (memory slots 1 to 30). A further press will show the overall average value out of all stored values. If a memory slot is empty, the symbol "-" will appear in the display for this slot.

Deleting stored data

The **medisana** Infrared clinical thermometer stores all measured values automatically. In case all 30 memory slots are being used, the oldest stored value will automatically be deleted to save a new result. If you want to delete all stored values, press and hold the **MEM**-button for at least 3 seconds, when the unit is switched on. The display will show the slot number "0" and "-". 4 beeps will be emitted for confirmation.

Error messages

Display	Description
Hi	The measured temperature is too high. Hi appears in the display if the reading is above 43°C (109.4°F) [body mode] or 100°C (212°F) [object mode].
Lo	The measured temperature is too low. Lo appears in the display if the reading is below 34°C (93.2°F) [body mode] or 0°C (32°F) [object mode].
Err	The ambient temperature is too high or too low. The ambient temperature is above 40°C (104°F) or below 16°C (60.8°F).
	Blank display. Check the position of the batteries.

Cleaning and maintenance

- The sensor is the most delicate part of the **medisana** infrared clinical thermometer. Protect it from dirt and damage. Clean the sensor with a cotton bud dipped in warm water or 70% isopropyl alcohol.
- Do not immerse the thermometer in water or any other liquid. **It is not watertight.**
- To clean the thermometer, use a soft cloth slightly dampened with mild soapy water or 70% isopropyl alcohol. Do not let water get into the unit. Only use the unit again when completely dry.
- Never use aggressive cleaning agents, thinners, benzene or tough brushes.
- Take care not to scratch the lens of the sensor, the display or the housing.
- Do not use or store the unit in an environment where the ambient temperature or humidity is too high or too low (see operating conditions in the "Technical specifications" section), in direct sunshine, in the vicinity of static electricity or in a dusty environment.
- Remove the batteries if you do not wish to use the unit for an extended period of time. Otherwise there is a danger that the batteries may leak.
- If the thermometer is used as described in these instructions, no regular maintenance or calibration is required.

Disposal

This product must not be disposed of together with domestic waste. All users are obliged to hand in all electrical or electronic devices, regardless of whether or not they contain toxic substances, at a municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner. Please remove the batteries before disposing of the device/unit. Do not dispose of old batteries with your household waste, but at a battery collection station at a recycling site or in a shop. Consult your municipal authority or your dealer for information about disposal.

Technical specifications

Name and model: **medisana** Infrared clinical thermometer **FTE**
digital display (LCD)
3 V =, 2 batteries 1.5 V (type LR03, AAA)
Body mode 34°C - 43°C (93.2°F - 109.4°F)
Object mode 0°C - 100°C (32.0°F - 212.0°F)
Body mode: In the region of 34°C - 35.9°C (93.2°F - 96.6°F): $\pm 0.3^{\circ}\text{C}$ (0.5°F)
36°C - 39°C (96.8°F - 102.2°F): $\pm 0.2^{\circ}\text{C}$ (0.4°F)
39.1°C - 43°C (102.4°F - 109.4°F): $\pm 0.3^{\circ}\text{C}$ (0.5°F)
(under operating conditions)
Object mode: $\leq 20^{\circ}\text{C}$ (68°F): $\pm 1^{\circ}\text{C}$ (1.8°F); $> 20^{\circ}\text{C}$ (68°F): $\pm 5\%$ 0.1°C (0.2°F)
up to 30 readings for temperatures $\geq 38^{\circ}\text{C}$ (100.4°F)
approx. 30 seconds after the last measurement 16°C - 40°C (60.8°F - 104°F), rel. humidity max. 85%
Storage/transport temperature: -20°C - +50°C (-4°F - +122°F), rel. humidity max. 85%
Dimensions: approx. 123 x 50.5 x 49 mm
Weight: approx. 80 g including batteries
Article No.: 48610
EAN Code: 40 15588 48610 4
CE 0297

In accordance with our policy of continual product improvement, we reserve the right to make technical and optical changes without notice.

The CE-mark of this device refers to the EC-directive 93/42/EEC MDD.

The unit is designed for home use only. For commercial use, a metrological inspection must be carried out once every year. This calibration will be charged for and can be carried out by an appropriate authority or an authorised service centre - in accordance with the regulations for users of medical equipment.

Electromagnetic compatibility:

The device complies with the EN 60601-1-2 standard for electromagnetic compatibility. Please check the appendix for details on this measurement data.

Warranty and repair terms

Please contact your dealer or the service centre in case of a claim under the warranty. If you have to return the unit, please enclose a copy of your receipt and state what the defect is.

The following warranty terms apply:

- The warranty period for **medisana** products is 3 years from date of purchase. In case of a warranty claim, the date of purchase has to be proven by means of the sales receipt or invoice.
- Defects in material or workmanship will be removed free of charge within the warranty period.
- Repairs under warranty do not extend the warranty period either for the unit or for the replacement parts.
- The following is excluded under the warranty:
 - All damage which has arisen due to improper treatment, e.g. nonobservance of the user instructions.
 - All damage which is due to repairs or tampering by the customer or unauthorised third parties.
 - Damage which has arisen during transport from the manufacturer to the consumer or during transport to the service centre.
 - Accessories which are subject to normal wear and tear.
- Liability for direct or indirect consequential losses caused by the unit are excluded even if the damage to the unit is accepted as a warranty claim.

medisana GmbH
Carl-Schurz-Str. 2, 41460 NEUSS, GERMANY

Electromagnetic compatibility - Guidance and manufacturer's declaration

Electromagnetic emissions		
The Infrared clinical thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of the Infrared clinical thermometer should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment-guidance
RF emissions CISPR 11	Group 1	The Infrared clinical thermometer uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Infrared clinical thermometer is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable	

Electromagnetic immunity

The **Infrared clinical thermometer** is intended for use in the electromagnetic environment specified below. The customer or the user of the **Infrared clinical thermometer** should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Electromagnetic immunity

The **Infrared clinical thermometer** is intended for use in the electromagnetic environment specified below. The customer or the user of the **Infrared clinical thermometer** should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the thermometer, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the thermometer is used exceeds the applicable RF compliance level above, the thermometer should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the thermometer.

b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the infrared clinical thermometer


The **Infrared clinical thermometer** is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the thermometer can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the **Infrared clinical thermometer** as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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