
## Unit and controls

- 1 Handle with electrodes
- 2 Trend light
- 3 VA LCD display
- 4 7-day trend graph
- 5 ITO tempered glass

### Legend:

 This instruction manual belongs to this device. The instruction manual includes information on the initial start-up and handling. Read this instruction manual completely. Failure to follow these in-structions may result in serious injury or damage to the equipment.

 **WARNING!**  
These warnings must be followed to pre-vent possible injury to the user.

 **NOTE**  
These notes give you useful additional information on the installation or operation.

 LOT number

 Manufacturer

# GB Safety instructions

 **Read the instruction manual carefully before using this device, especially the safety instructions, and keep the instruction manual for future use. Should you give this device to another person, it is vital that you also pass on these instructions for use.**

### INTENDED USE:

- This product is a device for estimating body composition by measuring the local bio-impedance value of the user's body. Used to measure body weight, body fat, body water, muscle mass, bone mass, BMI and basal metabolism, wireless transmission.
- This scale is not suitable for people with pacemakers or other medical implants. The indication of body fat may be inaccurate in the event of diabetes or other medical/physical conditions. The same applies to highly trained sportsmen.
  - The body fat measurement algorithm is calibrated for individuals between the ages of 6 and 80 years old. Users outside of this age range can only use the weight measurement feature.
  - The device uses Bioelectrical Impedance Analysis (BIA) and it may cause a weak current while taking measurement. So it's not recommended for pregnant women use.
  - Any treatment or diet for overweight or underweight requires qualified advice from appropriate specialists (doctor, dietician). The values determind with the scale can be of help.
  - Risk of tilting! Never place the scales on an uneven surface. Never step onto just one corner of the scale. While weighing, never stand on one side or on the edge of the scale.
  - Danger of slipping! Never step onto the scale with wet feet. Never step onto the scale while wearing socks.
  - Use the device only as directed by the instruction manual. In the event of misuse, all warrantees become null and void.
  - The scale is manufactured for home or personal use. It is not intended for commercial use in hospitals or other medical facilities.
  - Do not use the device if it does not work properly, if it has been dropped, dropped into water or damaged.
  - The scale has a measuring range up to 180 kg/396 lbs. Do not overload the scale.
  - Place the scale on a firm level surface. Soft, uneven surfaces are unsuitable for measurement and lead to wrong results.
  - Place the scale in a location where neither extreme temperatures nor extreme humidity are present.
  - Please handle the scale with care. Avoid shocks and vibrations of the scale. Do not drop it.
  - Stand carefully on the scale. Do not hop or jump on the treated surface, as this may damage the weighing mechanism.
  - Do not attempt to disassemble the device. Otherwise the warranty will become void. This device does not contain any parts that can be serviced or replaced by the user.
  - In the event of malfunctions, do not proceed with any unauthorised repair of the device. Otherwise, this will void any warranty claims. Have repairs performed only by an authorised service centre.

### Battery safety instructions

- Do not take the batteries apart!
- If required, clean the battery and device contacts before insertion!
- Removed spent batteries from the device immediately!
- Increased risk of leakage, avoid contact with skin, eyes and mucous membranes! In the event of contact with battery acid, immediately rinse the affected area with plenty of clean water and consult a doctor immediately!
- If a battery is swallowed, consult a doctor immediately!
- Insert the batteries correctly, paying attention to the polarity!
- Keep the battery compartment closed properly!
- Keep the batteries out of the reach of children!
- Do not re-charge the batteries! There is a risk of explosion!
- Do not short-circuit the device! There is a risk of explosion!
- Do not throw into a fire! There is a risk of explosion!
- Store unused batteries in the packaging and not near metal objects to avoid a short circuit!

### Scope of delivery

First, check whether the device is complete and is not damaged.  
If in doubt, do not operate the device and contact your service center.

- The scope of delivery includes:
- 1 **medisana** 56-in-1 Smart 8-electrode Body Analysis Scale
  - 1 Instruction manual
  - 4 x 1.5 V batteries (type AAA)


 Packaging is recyclable or can be recycled into raw materials. Please dispose of unwanted packaging material properly. If you notice any transport damage while unpacking, please contact your dealer immediately.

## WARNING

- Please ensure that the polythene packing is kept away from the reach of children!**
- Risk of suffocation!**
- The percentage of body fat, body water, muscle mass, bone weight, as well as body mass index (BMI), the values of the basal metabolic rate (KCAL) are to be considered only as approximate data. For any further information, refer to a doctor or to a dietician.**
- Self-measuring means control, it does not mean diagnosis or treatment. Any unusual levels should always be discussed with your doctor. Under no circumstances must you change the dose of medicines prescribed by your doctor.**

### How a body analysis scale works

**medisana** body analysis scale is tested by the Bioelectrical Impedance Analysis (BIA) method. The principle of Bioelectrical Impedance Analysis (BIA) is to apply weak current to body through the electrodes on the surface of scale, and use the theory that water in the body conducts electricity while fat does not, thereby measure the body's resistance. Then based on sample data, the body fat content is calculated by comparison.

**Low-frequency current:** obtain extracellular intelligence; high-frequency current: gain extracellular and Intracellular intelligence

**Fatty tissue:** Poor electrical conductivity

**Non-fatty tissue:** Better conductivity with more electrolyte

The body is divided into fatty tissue and non-fatty tissue according to above principle, and the resistance feedback of the two is also different. The body analysis scale calculates the proportion of different components of the body through data feedback and calculation.

### Requirements for correct measurement results

- Always weigh yourself in the same conditions.
- Measure at the same time of day.
- Stand upright and calm.
- It is best to take the measurement in the morning, after showering or bathing, with dry foot.

With the results of the measurement, when carried out consistently, a reliable control of the development of the body weight is possible.

### Insert/change batteries

Before putting your digital scale into operation, insert the four 1.5 V batteries (type AAA) included in the device. Open the battery compartment on the underside of the device and insert the batteries. Pay attention to the polarity (illustration in the battery compartment). Replace the battery compartment cover and push it in so that it clicks into place. Change the batteries if the battery change symbol "LO" appears on the display or if nothing is shown on the display after the unit has been operated.

### Initialise scale

Make sure that the scale is on firm and level ground . If you want to weigh yourself, but have moved the scale before, you must first initialise the scale. Briefly press the middle of the scale surface with your bare feet. The display will show "0.00". When the scale has turned off, it is ready for the "step-on" function. If the scale has not been moved before-hand, this process is omitted.

### Connect with the VigorHub App

#### View up to 56 body composition metrics via VigorHub App

- Download the VigorHub App free of charge from the App Store (iOS) or Google Play Store (Android).
- Enable Bluetooth and GPS on your mobile device. Open the VigorHub App and register a new account, then log in.
- From the dashboard, go to [Profile] > [Device Management] > [Add Device]. Follow the on-screen instructions to pair and connect your device.
- After connecting Bluetooth, press "Yes" to continue the Wi-Fi setup.
- Choose the Wi-Fi network and input password.
- Put your mobile device next to the scale. It will take around 5 seconds to complete the setup.

### User profile setup

For the first-time use, please set the personal information for the current user or new member in the VigorHub App.

### Change the weight unit

Tap profile and select [General settings] to change the weight unit directly in the VigorHub App.

### Body weighing and analysis

-  Stepping onto your monitor with socks or footwear on will result in an inaccurate measurement.
- If the posture is wrong, such as squatting, knees bent, the bare feet are not on the electrodes, the measurement will be inaccurate.

- Step onto the scale with bare feet, wrap your hands around the electrodes on both sides of the handle , and pull up the handle until your arms are level with your chest.
- Stand still and the final weight flashes on the display.
- Keep holding the handle and your body metrics will be gathered in around 3 seconds.
- When you see "-----" appears, this is to indicate that the data is being analyzed. Please maintain your posture during this period. Once the process is completed, the body weight, heart rate, body fat, body weight difference, muscle rate, bone mass, body water and BMI are displayed.
- The trend light  will remain blue during the weighing process. Once the weight is locked and user identification is complete, the light color will be fixed (red indicates heavier than last time, green indicates lighter or the same).
- The 7-day trend graph  will display the weight trend for the most recent 7 days based on the average value for each day.
- The scale can store 100 data points for 10 users on the scale. After the first measurement, the scale will be able to automatically recognize you based on the previous weight. If you measure offline, you will still see the 8 data points of your results. The data will be synchronized again once you reconnect to the app.

### Use Baby Holding mode via VigorHub App

You can measure baby's weight while holding the baby.

- Tap the " " icon on the dashboard. Tap [Device management] and select your paired device.
- Tap [Measure while holding baby]. Follow the guide to measure baby's weight.

### Activate Safe Pregnancy mode via VigorHub App

This body analysis scale calculates body indicators by measuring the safe current flowing through the body. It is not suitable for pregnant women and users with pacemakers.

Users can turn on/off the fat measurement on the App. Please follow the steps below:

- Tap the " " icon on the dashboard. Tap [Device management] and select your paired device.
- Tap the ":" icon at the top right corner and select [Switch mode].
- Turn on / off the fat measurement.

### Measurement symbols

|                                                                                       |                        |                                                                                                                                                                                                                                         |
|---------------------------------------------------------------------------------------|------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | Body weight            | Body weight refers to "the sum of all organs and metabolites in the body", which includes bones, muscles, fat, blood, internal organs, metabolic waste, and many other indicators. However, the largest proportion among them is water. |
|  | Heart rate             | Heart rate refers to the number of beats per minute of the heart. The normal heart rate range varies depending on individual differences.                                                                                               |
|  | Body fat               | Body fat percentage refers to the ratio of body weight in fatty tissue.<br>$\text{Body fat (\%)} = \frac{\text{Body fat weight (kg)}}{\text{Body weight (kg)}} \times 100\%$                                                            |
|  | Body weight difference | This shows the difference in body weight compared to the last measurement.                                                                                                                                                              |

|                                                                                     |             |                                                                                                                                                                                                                                                    |
|-------------------------------------------------------------------------------------|-------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | Muscle rate | Muscle rate refers to the percentage of muscle in the body composition to body weight. The higher the muscle rate, the greater the BMR, and the more calories consumed, making it less likely to gain weight.                                      |
|  | Bone mass   | Body bone mass refers to the total amount of bone tissue in the human Skeleton, including minerals and bone cells in bones. The size of bone mass directly affects bone health and plays a crucial role in the stability and function of the body. |
|  | Body water  | Body water ratio accounts for the largest ratio of body weight.<br>$\text{Body water (\%)} = \frac{\text{Weight of body water}}{\text{Weight of body fluid}} \times 100\%$                                                                         |
|  | BMI         | Body Mass Index (BMI) is the relationship between weight and height, body fat and health risk.<br>$\text{BMI} = \frac{\text{Body weight (kg)}}{\text{Body height (m)}^2}$                                                                          |

### Error messages

| Abnormal displays | Cause                                      | Solution                                                                                       |
|-------------------|--------------------------------------------|------------------------------------------------------------------------------------------------|
| Err0              | Bioelectrical impedance measurement failed | Please check whether the scale is put freely on a firm and level surface.                      |
| Err1              | Resting heart rate measurement failed      | It must not touch a wall or another object. Keep your body still with bare feet and try again. |
| Err2              | Failed to turn on Wi-Fi                    | Please check the network connection and try again.                                             |
| Err3              | Router connection failed                   |                                                                                                |
| Err4              | Wi-Fi upload failed                        |                                                                                                |
| Err5              | Wi-Fi distribution network failed          |                                                                                                |

### Cleaning and care

Never use aggressive cleaning agents or stiff brushes. Clean the scale only with a soft, slightly damp cloth. Never use harsh detergents or alco-hol. No water should enter the device. Only use the device again when it is completely dry.

### Instructions for disposal

-  This device must not be disposed of with domestic waste. Metals are reusable or can be recycled into raw materials. Dispose of a device that is no longer needed at a collection point in your town/city so that it can be disposed of in an environmentally friendly way. Please contact your local authority or retailer with regard to disposal procedures. Remove the battery before you dispose of the device. Do not dispose of spent batteries in domestic waste, rather in special waste or at a battery collection point at your retailer!

### Technical specifications

|                                   |                                                                                           |
|-----------------------------------|-------------------------------------------------------------------------------------------|
| Name and model no.                | : <b>medisana</b> 56-in-1 Smart 8-electrode Body Analysis Scale<br><b>BS 860 connect</b>  |
| Power supply                      | : 6 V (4 x 1.5V, AAA)                                                                     |
| Units of measurement              | : kg / st / lb                                                                            |
| Weighing range                    | : 0.1 - 220 kg                                                                            |
| Division                          | : 0.05 kg (<100 kg); 0.1 kg (>100 kg)                                                     |
| Weight measurement accuracy       | : ±0.1 kg (0.1-50 kg); ±0.2 kg (50-100 kg);<br>±0.3 kg (100-150 kg); ±0.5 kg (150-220 kg) |
| Operating conditions              | : +5°C to +40°C, ≤ 90% relative humidity                                                  |
| Storage/Transportation conditions | : -20°C to +50°C, ≤ 90% relative humidity                                                 |
| Dimensions                        | : 350 x 350 x 25 mm                                                                       |
| Weight                            | : 2.4 kg                                                                                  |
| Article no.                       | : 50059                                                                                   |
| EAN no.                           | : 4897138500597                                                                           |

**In the course of constant product improvements, we reserve the right to make technical and design changes.**

### Warranty and repair conditions

In case of warranty please contact your specialist shop or the service centre directly. If you need to return the device, please indicate the defect and enclose a copy the purchase receipt. The following warranty conditions apply:

- All **medisana** products are guaranteed for one year from the date of purchase.  
The date of purchase is to be proven in case of warranty by the purchase receipt or invoice.
- Defects due to material or manufacturing defects shall be eliminated free of charge within the warranty period.
- A warranty service does not extend the warranty period for either the equipment or any replacement parts.
- The following are excluded from the guarantee:
  - any damage caused by improper handling, e.g. by non-observance of the instruction manual.
  - Damage due to repair or intervention by the purchaser or unauthorised third parties.
  - Transport damage that has occurred on the way from the manufacturer to the consumer or when sending it to the service centre.
- Liability for direct or indirect consequential damages caused by the device is also excluded if the damage to the device is recognised as a warranty case.

|                                                                                    |                        |                                                                                     |                        |                                                                                       |                       |                                                                                       |                               |
|------------------------------------------------------------------------------------|------------------------|-------------------------------------------------------------------------------------|------------------------|---------------------------------------------------------------------------------------|-----------------------|---------------------------------------------------------------------------------------|-------------------------------|
|    | Weight                 |    | Left arm fat ratio     |    | Left leg muscle mass  |    | Weight control                |
|    | Heart rate             |    | Right arm fat ratio    |    | Right leg muscle mass |    | Standard weight               |
|    | BMR                    |    | Left arm fat mass      |    | Trunk fat percentage  |    | Fat-free mass                 |
|    | BMI                    |    | Right arm fat mass     |    | Trunk fat mass        |    | Visceral fat                  |
|    | Fat ratio              |    | Left arm muscle rate   |    | Trunk muscle rate     |    | Bone mass                     |
|    | Fat mass               |    | Right arm muscle rate  |    | Trunk muscle mass     |    | Skeletal muscle quality index |
|   | Muscle rate            |   | Left arm muscle mass   |   | Intracellular water   |   | Body type                     |
|  | Muscle mass            |  | Right arm muscle mass  |  | Extracellular water   |  | Body age                      |
|  | Protein ratio          |  | Left leg fat ratio     |  | Skeletal muscle ratio |  | Health evaluation             |
|  | Protein mass           |  | Right leg fat ratio    |  | Skeletal muscle mass  |  | Health score                  |
|  | Water ratio            |  | Left leg fat mass      |  | Obesity level         |  | Body cell mass                |
|  | Total body water       |  | Right leg fat mass     |  | Adiposity level       |  | Minerals                      |
|  | Subcutaneous fat ratio |  | Left foot muscle rate  |  | Fat control           |  | Waist-hip ratio               |
|  | Subcutaneous fat mass  |  | Right foot muscle rate |  | Muscle control        |  | Recommended calorie intake    |

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